
ORTHOPAEDICS

Instructions for Use

Mitchell Ponseti® Ankle Foot Orthosis (AFO) and Ponseti® Abduction Bar

MD Orthopaedics, Inc.
604 North Parkway Street
Wayland, IA 52654 USA
www.mdorthopaedics.com

Document MD-124
Version 8.0
January 2022

Emergo Europe B.V.
Prinsessegracht 20, 2514 AP
The Hague, The Netherlands

MedEnvoy
Prinses Margrietplantsoen 33
Suite 123
2595 AM The Hague
The Netherlands

EN

Table of Contents

| | |
|----------------------------------------------------------------------------------------------|-----------|
| DISCLAIMER..... | 1 |
| TRADEMARKS | 1 |
| 1 PREFACE | 1 |
| 1.1 Intended Purpose | 1 |
| 1.2 Target Patient Group and Intended Users | 2 |
| 1.3 Clinical Benefits | 2 |
| 1.4 Explanation of Safety Warnings | 2 |
| 1.5 Retention Instructions..... | 2 |
| 1.6 Obtaining Documentation and Information | 2 |
| 1.6.1 Ordering Documentation..... | 2 |
| 1.6.2 Other languages | 2 |
| 1.6.3 Documentation Feedback..... | 2 |
| 1.6.4 Support and service..... | 2 |
| 1.6.5 Name and address of the manufacturer | 3 |
| 2 DESCRIPTION OF THE PRODUCT | 3 |
| 2.1 Intended Use and Reasonably Foreseeable Misuse..... | 3 |
| 2.2 Sterilization State and Method | 3 |
| 2.3 Summary of Safety and Clinical Performance | 3 |
| 2.4 Technical Characteristics..... | 4 |
| 2.4.1 Mitchell Ponseti® Ankle Foot Orthotic..... | 5 |
| 2.4.2 Ponseti® Abduction Bar..... | 5 |
| 3 SAFETY INSTRUCTIONS | 5 |
| 3.1 How to Use the Product Safely | 6 |
| 3.1.1 Technical life span and Warranty..... | 6 |
| 3.1.2 Safety information related to the intended use and reasonably foreseeable misuse | 6 |
| 3.1.3 Product limitations and restrictions and contraindications | 6 |
| 3.1.4 Safety information when using the device in combination with other devices..... | 6 |
| 3.1.5 Safety information regarding the use | 6 |
| 3.1.6 Safe Disposal | 6 |
| 3.2 Potential Health Consequences..... | 6 |
| 4 INSTRUCTIONS FOR USE..... | 7 |
| 4.1 Clubfoot Treatment Overview | 7 |
| 4.2 Mitchell Ponseti® AFO..... | 8 |
| 4.3 Ponseti® Abduction Bar | 9 |
| 5 PREPARATION | 9 |
| 5.1 How to Transport and Store the Product..... | 9 |
| 6 MAINTENANCE | 9 |
| 6.1 Reusing the Device | 9 |
| 6.1.1 Cleaning the device..... | 9 |
| 6.2 How to Inspect the Product | 9 |
| 7 TROUBLESHOOTING | 10 |
| 7.1 How to Identify and Solve Problems | 10 |
| 7.1.1 Troubleshooting by non-skilled persons..... | 10 |
| 7.2 Frequently Asked Questions..... | 10 |
| 8 GLOSSARY | 11 |
| 9 ICON LEGEND | 12 |

DISCLAIMER

MD Orthopaedics makes no representations or warranties with respect to this manual and, to the maximum extent permitted by law, expressly limits its liability for breach of any warranty that may be implied to the replacement of this manual with another. Furthermore, MD Orthopaedics reserves the right to revise this publication at any time without incurring an obligation to notify any person of the revision.

As the designer and manufacturer of products, MD Orthopaedics does not provide medical treatment and/or advice. Information about the risks and benefits of medical treatment, including treatment that involves the use of MD Orthopaedics products and available alternatives is provided solely by physicians or other health care providers. If you have questions about your child's treatment, it is important to discuss those questions with the appropriate health care provider.

MD Orthopaedics does not, and cannot, warrant certain results from the use of its products.

The information provided in this documentation contains general descriptions and/or technical characteristics of the performance of the products contained herein. This documentation is not intended as a substitute for and is not to be used for determining suitability or reliability of these products for specific user applications. It is the duty of any such user to perform the appropriate and complete risk analysis, evaluation, and testing of the products with respect to the relevant specific application or use thereof. Neither MD Orthopaedics nor any of its affiliates or subsidiaries shall be responsible or liable for misuse of the information that is contained herein. If you have any suggestions for improvements or amendments or have found errors in this publication, please notify us.

All pertinent state, regional, and local safety regulations must be observed when using this product. For reasons of safety and to help ensure compliance with documented system data, only the manufacturer shall perform repairs to components.

When devices are used for applications with technical safety requirements, the relevant instructions must be followed. Failure to observe this information can result in injury or equipment damage.

MD Orthopaedics has made every effort to prevent that the labelling, instructions for use, making available, putting into service and advertising of devices uses text, names, trademarks, pictures and figurative or other signs that may mislead the user or the patient regarding the device's intended purpose, safety, and performance.

Copyright © 2022 by MD Orthopaedics, Inc.

All rights reserved. No part of this publication may be reproduced, distributed, or transmitted in any form or by any means, including photocopying, recording, or other electronic or mechanical methods, without the prior written permission of the publisher. For permission requests, write to the publisher at the address below.

MD Orthopaedics, Inc.
604 North Parkway Street
Wayland, IA 52654
USA
+1-877-766-7384
www.mdorthopaedics.com

TRADEMARKS

Mitchell Ponseti® and **Ponseti**® are registered trademarks of MD Orthopaedics.

MD Orthopaedics has made every effort to supply trademark information about company names, products and services mentioned in this manual. Trademarks shown below were derived from various sources. All trademarks are the property of their respective owners.

General Notice: Some product names used in this manual are used for identification purposes only and may be trademarks of their respective companies.

1 PREFACE

1.1 Intended Purpose

The **Mitchell Ponseti**® **Ankle Foot Orthosis (AFO)** and **Ponseti**® **Abduction Bar** is an ankle-foot orthosis, also known as **Mitchell Ponseti**® **Brace**, intended for use in the Ponseti Method of treatment for the correction of Congenital Talipes Equinovarus (CTEV) in infants and children less than or equal to eight years of age. The brace is intended to prevent relapse (equinus and varus deformity of the heel).

The **Mitchell Ponseti**® **Brace** is to be prescribed by a physician or health care provider trained in the Ponseti Method of clubfoot treatment and is typically worn for 23 hours a day for three months after casting and then during night-time and nap time for 4-5 years or until the physician or health care provider indicate the device is no longer needed.

A clinical orthotist, also trained in the Ponseti Method, may assist the physician or health care provider in determining the correct size of the AFO for the child. The orthotist may also train the child's parents and families in the right way to put the AFO onto the child's foot, what cues may indicate poor fit or other issues, and how to determine it is time to move to the next size.

Other members of the clinical team trained in the Ponseti Method – nurses, physician assistants, or other licensed healthcare professionals – may also train parents on the proper use of the **Mitchell Ponseti**® **Brace** and make any necessary adjustments in fit or sizing. The child's parents or family or other designated caretaker are responsible for adhering to the prescribed Ponseti Method treatment regimen of 23 hours a day for three months, then during sleep and during naptime for the next 4-5 years or until the physician indicates the device is no longer needed.

1.2 Target Patient Group and Intended Users

The MD Orthopaedics **Mitchell Ponseti® Ankle Foot Orthosis (AFO)** and **Ponseti® Abduction Bar** is designed to be worn by a child (≤8) years of age with clubfoot after the casting portion of treatment is complete.

This document is intended for the physician, orthotist, clinical team, parents, and families of children receiving treatment with the **Mitchell Ponseti® Ankle Foot Orthosis (AFO)** and **Ponseti® Abduction Bar**.

1.3 Clinical Benefits

The Ponseti Method is a non-surgical clubfoot treatment process developed by Dr. Ignacio Ponseti. The method utilizes gentle manipulation and casting to gradually move the foot into the correct position. Each week the cast is removed, the foot carefully manipulated into the next position, and a new cast is applied. This process continues for up to 6 to 8 weeks.

Following the casting period, the feet are held in the proper position using a foot abduction brace (FAB). The brace is comprised of the Ankle Foot Orthosis (AFO) and the connecting, abduction bar. The brace is typically worn for 23 hours a day for three months, then worn at night and during naptime for up to 4-5 years or as prescribed by a physician to ensure the correction holds.

With proper treatment and adherence to the bracing protocol prescribed by the physician, it is expected that most cases of clubfoot can be corrected by the time the child reaches 4 or 5 years of age.

MD Orthopaedics does not provide medical treatment and/or advice. Information about the risks and benefits of medical treatment, including treatment that involves the use of MD Orthopaedics products should be provided solely by the physician or other qualified health care providers treating your child. If you have questions about your child's treatment, it is important for you to discuss those questions with the appropriate health care provider.

1.4 Explanation of Safety Warnings

Caution indicates a hazard with a low level of risk which, if not avoided, could result in minor or moderate injury.

Indicates information considered important, but not hazard-related.

1.5 Retention Instructions

Read and understand this manual and its safety instructions before using this product. Failure to do so can result in injury.

Follow all the instructions.

The product shall only be used by persons who have fully read and understand the contents of this user manual.

Ensure that each person who uses the product has read these warnings and instructions and follows them.

The manufacturer is not liable for cases of material damage or personal injury caused by incorrect handling or non-compliance with the safety instructions. In such cases, the warranty will be voided.

This policy / warranty is non-transferable and covers normal wear and tear only. Altering the product and after-market modifications will void the warranty. Please contact us if you experience problems not outlined in this policy. Atypical matters will be handled on a case-by-case basis.

1.6 Obtaining Documentation and Information

1.6.1 Ordering Documentation

Additional documentation, user instructions and technical information can be ordered by calling MD Orthopaedics at +1-877-766-7384 or e-mailing info@mdorthopaedics.com.

1.6.2 Other languages

Instructions for Use (IFU) manuals are available in other languages at www.mdorthopaedicsifu.com.

1.6.3 Documentation Feedback

If you are reading MD Orthopaedics' product documentation on the internet, any comments can be submitted on the support website. Comments can also be sent to info@mdorthopaedics.com. We appreciate your comments.

1.6.4 Support and service

For questions, information, technical assistance or to order user instructions, please contact:

MD Orthopaedics, Inc.

604 North Parkway Street

Wayland, IA 52654

USA

www.mdorthopaedics.com

+1-877-766-7384 or e-mail info@mdorthopaedics.com

1.6.5 Name and address of the manufacturer

The following natural or legal person makes the device, to which this user manual applies, suitable for use within the European Union and is the manufacturer of the device:

The **Mitchell Ponseti® Ankle Foot Orthosis** and **Ponseti® Abduction Bar** are manufactured by:

MD Orthopaedics, Inc.
604 North Parkway Street
Wayland, IA 52654
USA
www.mdorthopaedics.com
+1-877-766-7384

2 DESCRIPTION OF THE PRODUCT

2.1 Intended Use and Reasonably Foreseeable Misuse

The **Mitchell Ponseti® Brace** is a Class I medical device intended for use only on children whose clubfoot is being corrected with the Ponseti Method of treatment, and only as prescribed by a physician fully trained in the Ponseti Method.

The **Mitchell Ponseti® Brace** shall not be used without prescription by a physician trained in the Ponseti method and is to be used only as prescribed.

The **Mitchell Ponseti® Brace** shall not be used to treat any other orthopedic condition other than clubfoot unless otherwise prescribed by the physician or orthotist.

The **Mitchell Ponseti® Ankle Foot Orthosis** shall be used with the following original accessories and components only:

- **Ponseti® Adjustable Bar**
 - Available with either 10 or 15 degrees of dorsiflexion.
- **Pressure Saddle:** This soft, comfortable pressure saddle attaches to the middle strap of the **Mitchell Ponseti® AFO** to provide additional relief of pressure areas. Available in size small and large. Size small is recommended for use with AFO Size 0000-1; Size large is recommended for use with Size 2-12.
- **Bar Cover:** This soft bar cover provides protection and comfort for your child, and is available in size small (5") and large (8"), in pink, light blue, or dark blue.

The **Mitchell Ponseti® Ankle Foot Orthosis** is available in four styles:

The **Mitchell Ponseti® Standard AFO** uses the Ponseti Method of clubfoot correction. Our patented system features a soft lining for comfort and compliance, soft synthetic leather straps. The location of the heel can easily be seen through two viewing holes in the back of the AFO. The AFO is worn attached to the **Ponseti® Abduction Bar** in accordance with the standard Ponseti bracing protocol. Available in sizes preemie through 12.

The **Mitchell Ponseti® Plantar Flexion Stop AFO:** The articulated Plantar Flexion Stop limits plantar flexion with rigid dorsiflexion to maintain the prescribed position of the foot. Recommended for the complex clubfoot. This device is also good for children with hyper flexible feet as it ensures the foot is always held in a neutral or dorsiflexed position. The AFO is worn attached to the **Ponseti® Abduction Bar** in accordance with the standard Ponseti bracing protocol. Available in sizes 00 through 12.

The **Mitchell Ponseti® Toe Stilt AFO** was specifically designed by recommendation of Dr. Ponseti to maintain clinically desirable dorsiflexion by stretching the Achilles tendon and helping to keep the foot flexible. This is achieved by the placement of a built in wedge attached to the front of the AFO which keeps both feet in a dorsiflexed position of 10 degrees when standing. The AFO is worn attached to the **Ponseti® Abduction Bar** in accordance with the standard Ponseti bracing protocol. The patient must be able to stand independently for the Toe Stilt AFO to have any effect when used without the **Ponseti® Abduction Bar**. Available in sizes 2 through 11.

The **Mitchell Ponseti® Plantar Flexion Stop/Toe Stilt** is a dual purpose AFO that supports the Ponseti Method of treating clubfoot. The AFO is worn attached to the **Ponseti® Abduction Bar** in accordance with the standard Ponseti bracing protocol. The Plantar Flexion Stop (PFS) is designed to structurally limit any plantar flexion and work particularly well in clubfoot cases with hypermobile or corrected atypical/complex deformity. The PFS is mounted on the standard AFO during the assembly process but is designed to be custom fit as necessary by the orthotist to meet the dynamic needs of each clubfoot child. The Toe Stilt encourages lower calf, Achilles tendon and plantar fascia stretch in the toddler who can stand independently without the **Ponseti® Abduction Bar**. The Toe Stilt holds both feet at 10 degrees dorsiflexion when standing. Available in sizes 2 through 11.

 Education on proper brace application and troubleshooting for skin issues is critical to ensure your child's clubfoot correction is effective. DO NOT USE this device without first being trained by your physician or clinical team. Contact Customer Support for assistance in locating training resources.

2.2 Sterilization State and Method

Not applicable, as the **Mitchell Ponseti® Brace** is not a sterile product.

2.3 Summary of Safety and Clinical Performance

The **Mitchell Ponseti® Brace** is a Class I, non-implantable device. Therefore, a Summary of Safety and Clinical Performance is not required.

2.4 Technical Characteristics

The **Mitchell Ponseti® AFO** was designed with a soft, contoured elastomeric liner to cushion and hold the child's foot in place. Straps and a body made from a ductile synthetic suede conform to the shape of the foot for added comfort and stability. The AFO has openings on the side to add ventilation and keep the foot cool, and viewing holes to help gauge proper placement of the heel in the AFO.

NOTE: Socks are **required** to eliminate skin contact with product.

The **Ponseti® Abduction Bar** was designed to be easily adjustable to the shoulder width of the child. The Quick Clips make it easy for the bar to be attached and removed from the AFOs as needed.

The length of the foot should be measured from the heel to tip of big toe. Do not add growing room as this is allowed for in the table below.

The Premie AFOs are permanently attached to an adjustable bar.


MD Orthopaedics' AFO Sizing Chart

| Mitchell Ponseti® AFO Size | (A) Foot Length (cm) | (B) Foot Width (cm) |
|----------------------------|----------------------|---------------------|
| P6-0 * | 5.5 to 6.5 | 2.5 to 4.0 |
| P5-0 * | 5.8 to 6.9 | 3.0 to 4.5 |
| 0000 | 6.0 to 7.3 | 3.3 to 4.8 |
| 000 | 6.3 to 7.8 | 3.0 to 5.5 |
| 00 | 6.8 to 8.3 | 3.0 to 5.5 |
| 0 | 7.5 to 9.0 | 3.0 to 6.0 |
| 1 | 8.0 to 9.5 | 3.5 to 6.0 |
| 2 | 9.0 to 10.5 | 3.8 to 6.5 |
| 3 | 10.0 to 11.5 | 4.5 to 7.0 |
| 4 | 11.0 to 12.5 | 4.7 to 7.0 |
| 5 | 12.0 to 13.5 | 4.8 to 7.0 |
| 6 | 13.0 to 14.5 | 5.0 to 7.5 |
| 7 | 14.0 to 15.5 | 5.5 to 8.0 |
| 8 | 15.0 to 16.5 | 6.7 to 8.2 |
| 9 | 16.0 to 17.5 | 7.0 to 8.5 |
| 10 | 17.0 to 18.5 | 7.4 to 8.9 |
| 11 | 18.0 to 19.5 | 7.7 to 9.2 |
| 12 | 19.0 to 20.5 | 8.0 to 9.5 |

Legacy European AFO Sizing Chart

| Mitchell Ponseti® AFO Size | (A) Foot Length (mm) |
|----------------------------|----------------------|
| P6-0 (Premie 1) | 60 to 62 |
| P5-0 (Premie 2) | 63 to 66 |
| 0000 | 67 to 70 |
| 000 | 71 to 75 |
| 00 | 76 to 80 |
| 0 | 81 to 87 |
| 1 | 88 to 92 |
| 2 | 93 to 102 |
| 3 | 103 to 112 |
| 4 | 113 to 122 |
| 5 | 123 to 132 |
| 6 | 133 to 142 |
| 7 | 143 to 152 |
| 8 | 153 to 162 |
| 9 | 163 to 172 |
| 10 | 173 to 182 |
| 11 | 183 to 192 |
| 12 | 193 to 202 |

Length and width measurement locations

| Ponseti® Abduction Bar Size ** | Minimum Width (cm) | Maximum Width (cm) |
|--------------------------------|--------------------|--------------------|
| Extra Short *** | 17.3 | 25.0 |
| Short | 20.0 | 30.0 |
| Long | 23.5 | 37.8 |


When sizing, we recommend leaving 0.5-1.0 cm room for growth; i.e. a foot length of 9.0 cm would be a size 1. Please email us for sizing assistance at info@mdorthopaedics.com.

*The **Mitchell Ponseti®** Size P6-0 or P5-0 AFOs are permanently attached to a Premie Bar, which is pre-set to 10° dorsiflexion and 60° abduction. The bar adjusts in width from 15.0-20.0 cm, in 1 cm increments. The preemie bar can be cut down as small as 10 cm, upon request.


The **Ponseti® Abduction Bars are available in 10 or 15° dorsiflexion.

***The **Ponseti® Abduction Bar**, size Extra Short Bar can be cut down as small as 11 cm, upon request.

2.4.1 Mitchell Ponseti® Ankle Foot Orthotic

2.4.2 Ponseti® Abduction Bar

A hex key (wrench) is also included, for making adjustments of the abduction bar as described in section 4.3.

3 SAFETY INSTRUCTIONS

Use the **Mitchell Ponseti® Ankle Foot Orthosis (AFO)** and **Ponseti® Abduction Bar** only in accordance with these instructions and as directed by your physician, health care provider or other member of the clinical team.

Improper application of the AFO or bar can lead to skin irritation, sores, and blisters.

Socks are **required** to eliminate skin contact with product and prevent skin irritation and blisters. Pressure Saddles are recommended.

Consult your physician, health care provider or clinical team immediately if the child is experiencing pain, blisters, sores, or if the AFO is visibly slipping on the child's foot.

Do not bend the bar or allow children to walk with the bar attached. Bending of the bar may compromise therapeutic value and lead to breakage in use, which may cause harm to the patient.

IMPORTANT: Ensure the dressing area is well lit and free from distraction to ensure you apply the AFO properly. Improper application of the AFO or bar can lead to skin irritation, sores, and blisters.

⚠️ Read and understand this manual and its safety instructions before using this product. Failure to do so can result in injury.

3.1 How to Use the Product Safely

3.1.1 Technical life span and Warranty

- On average, each AFO will last 3-9 months, depending on the child's growth rate. As the child's growth rate slows, more time can be expected.
- The bars are adjustable and may last more than one pair of AFOs if adjustable bar width continues to facilitate shoulder width requirement. The bar width should be periodically checked and adjusted to match the child's shoulder width or as recommended by your physician, orthotist or health care provider.
- Warranty of the **Mitchell Ponseti® Brace** is one year under normal use, and is single patient multiple use. Normal use is defined by Normal Daily Use for 365 days/year.

3.1.2 Safety information related to the intended use and reasonably foreseeable misuse

- Use the **Mitchell Ponseti® Brace** only in accordance with these instructions and as directed by your physician or other member of the clinical team.
- Consult your physician or clinical team immediately if the child is experiencing pain, blisters, sores, or if they are not in good condition or do not correctly fit the child.
- Do not overtighten the adjusting screws.
- Do not use the **Mitchell Ponseti® Brace** if it is not in good condition or does not correctly fit the child.
- Protect your child, yourself, and your furniture from being hit by the bar when the child is wearing it. It is recommended to pad the bar with a **Ponseti® Bar Cover**.

3.1.3 Product limitations and restrictions and contraindications

- Do not bend the bar or allow children to walk with the bar attached. Bending of the bar may compromise therapeutic value and lead to breakage in use, which may cause harm to the patient.
- Do not use the **Mitchell Ponseti® Brace** if it is not in good condition or does not correctly fit the child.
- The **Mitchell Ponseti® Brace** shall not be used without prescription by a physician or health care provider trained in the Ponseti Method and is to be used only as prescribed.
- The **Mitchell Ponseti® Brace** shall not be used to treat any other orthopedic condition other than clubfoot unless otherwise prescribed by the physician or health care provider.

3.1.4 Safety information when using the device in combination with other devices

- The **Mitchell Ponseti® Ankle Foot Orthosis** shall be used with the following original accessories and components only:
 - **Ponseti® Abduction Bar**
 - **Pressure Saddle**
 - **Ponseti® Bar Cover**

 Altering the product and/or making after-market modifications will void the warranty and may result in injury.

3.1.5 Safety information regarding the use

- Consult your physician or clinical team immediately if the child is experiencing pain, blisters, sores, or if the AFO is visibly slipping on the child's foot.
- Do not bend the bar or allow children to walk with the bar attached. Bending of the bar may compromise therapeutic value and lead to breakage in use, which may cause harm to the patient.

3.1.6 Safe Disposal

- The **Mitchell Ponseti® Brace** and accessories can be disposed of in the regular trash. None of the components of the device are recyclable.

3.2 Potential Health Consequences

It is necessary that the brace be worn 23 hours a day for three months and then 12-14 hours a day (naps and nighttime) for 4-5 years or as directed by your physician or health care provider.

Ensure the area is well lit and free from distraction to ensure you apply the AFO properly. Improper application of the AFO or bar can lead to skin irritation, sores, and blisters.

Incorrect use of the AFO may result in relapse of the clubfoot defect or injury to the child's foot. If a user and/or patient experiences any serious incident that has occurred in relation to the device, it should be reported to the MD Orthopaedics, Inc. and the competent authority of the Member State in which the user and/or patient is established.

Please contact MD Orthopaedics if you experience problems not outlined in this policy.


MD Orthopaedics, Inc.
604 North Parkway Street
Wayland, IA 52654
USA

www.mdorthopaedics.com

+1-877-766-7384 or e-mail info@mdorthopaedics.com

4 INSTRUCTIONS FOR USE

4.1 Clubfoot Treatment Overview

This infographic provides information — not advice. If you need medical advice, please consult a doctor or other appropriate medical professional.


No warranties are given in relationship to the medical information supplied and no liability will be accrued to MD Orthopaedics in the event the user suffers loss because of reliance upon the information provided in this infographic.

4.2 Mitchell Ponseti® AFO

Carefully inspect the AFO before each use. Never use an AFO or bar that is damaged, broken or not functioning properly, as these may injure your child and/or render treatment ineffective. Always use socks with the AFO to prevent direct contact with the skin.

1. Open the AFO so all the suede straps are clear for inserting the foot. Buckles are designed to be placed to the inside of the foot. They may be positioned to the outside based on physician preference or individual convenience of application.


2. Holding the lower leg, gently slide the foot into place until the heel is secure against the back and bottom of the AFO. Socks that cover the foot and lower leg are required to prevent skin contact with the AFO.


The exact position of the heel in the AFO will vary depending on the type and severity of the clubfoot. Forcing the heel into position can lead to sores, so check with your physician to verify the proper position.

3. Pull tongue horizontally against the ankle and hold in place with your thumb. Be sure the hole in the tongue is set in the center of the ankle, above the middle strap.

4. Buckle the middle strap securely over the horizontal tongue. Securely buckle the ankle strap but be careful not to overtighten, as that can lead to skin irritation and sores.

5. Look through the heel viewing holes at the back of the AFO to ensure the heel is ideally down and to the back of the AFO. The heel may not initially be able to contact the bottom of the AFO after the tenotomy and last casting but will over time. If you do not see the heel in the viewing holes, readjust the tongue and middle strap so that the heel is properly positioned.

In some cases, the heel will not touch the bottom of the AFO, and in others the heel will be barely visible through the heel viewing holes. In these challenging cases, the heel will drop over time as long as ankle motion is adequate. Forcing the heel into the AFO and overtightening the middle strap can lead to skin irritation, sores, and reduced bracing compliance.


6. Once upper straps are tight and heel is securely in place, buckle the toe strap and retighten the other straps if needed.


4.3 Ponseti® Abduction Bar

The bar should be set so that the width of the bar is equal to - or even slightly wider than - the shoulder width of the child. Measure the shoulder width of the child from the outside left shoulder to the right outside shoulder. Adjust the length of the bar, measuring from the center screw of the left heel cap to the center screw of the right heel cap, so that it is equal to the shoulder width measurement. It is important to measure the shoulder width regularly and adjustments should be approved by orthotist.

1. Use the enclosed hex key to loosen set screws in bar lock. Adjust the width of the bar to the child's shoulder width, using the center screw in the heel caps as a guide. Re-tighten the set screws in the bar lock.

2. Using a screwdriver, loosen center screw of heel cap; remove outer screw. Rotate clip to desired degree indicator to set abduction (external rotation) as determined by your physician. Re-insert outer screw into appropriate hole. Tighten both screws.

3. Insert clip into sole of AFOs until you hear a "click" - there will be an audible click when they are seated properly. Verify that the bar lock and Quick Clips are engaged when assembled by pulling on the AFO to make sure it does not disconnect from the bar. If damage is found discontinue use and contact Customer Service.

5 PREPARATION

5.1 How to Transport and Store the Product

The **Mitchell Ponseti® Ankle Foot Orthosis** and the **Ponseti® Abduction Bar** ship at room temperature and require no special handling during transport or storage.

6 MAINTENANCE

6.1 Reusing the Device

All parts were designed to be used by a single patient under normal use conditions and as prescribed by a physician.

Normal use includes the daily application of the **Mitchell Ponseti® Brace** and use up to 23 hours per day, or as prescribed.

6.1.1 Cleaning the device

To clean the device:

1. Hand or machine-wash the AFO using cool water, a gentle cycle and mild detergent.
2. Let the AFO air dry completely before use.
3. It is not necessary to wash before using.

 Do not use hot water when washing the **Mitchell Ponseti® Ankle Foot Orthosis**. This may cause the synthetic suede straps to warp, shrink, separate or delaminate. Do not put the AFO into a mechanical dryer. Air dry only.

6.2 How to Inspect the Product

The **Mitchell Ponseti® Ankle Foot Orthosis (AFO)** and **Ponseti® Abduction Bar** should be inspected every time for proper fit, manufacturing defect, and damage prior to the device being placed onto the child's feet.

 Carefully inspect the **Mitchell Ponseti® Ankle Foot Orthosis (AFO)** and **Ponseti® Abduction Bar** before each use. Never use an AFO or bar that is damaged, broken or not functioning properly, as this may injure your child and/or render treatment ineffective. Contact Customer Service for Support, if needed.

7 TROUBLESHOOTING

7.1 How to Identify and Solve Problems

7.1.1 Troubleshooting by non-skilled persons

| Error | Cause | Solution |
|-----------------------------------------------------|------------------------------------------------|-----------------------------------------------|
| Child develops blisters or other damage to the skin | Improper fit and/or use | Contact your physician or healthcare provider |
| Broken AFO (within 6 months of wear) | Manufacturing defect or damaged via normal use | Contact your point of purchase |
| Broken AFO (after 6 months of wear) | Damaged via normal use | Contact your point of purchase |
| Broken Bar, Quick Clip, Heel Cap or Bar Lock | Damaged via normal use | Contact your point of purchase |
| Incorrect Size (within 7 days of receipt) | Sizing Error | Contact your point of purchase |

7.2 Frequently Asked Questions

| Question | Answer |
|------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Do the buckles go to the inside of the foot or to the outside of the foot? | The AFO was designed with the buckles set to the inside of the foot. This design feature is to aid parents in the donning and doffing of the AFOs. Ponseti experts do allow buckles to the outside of the foot based on professional preference or customer convenience. It will not alter the performance or risk of injury by orienting the buckles outward. |
| I just received my order, and the AFOs seem too big/too small. Is there a different size that would work better? | Contact your point of purchase |
| The child's heel is not down in the AFO. Is something wrong? | If the heel is not visible in the viewing holes, check with your physician to ensure the clubfoot is fully corrected. If the physician indicates that the correction is good, it may not be necessary to force the heel down into the AFO, as this can lead to skin irritation, bracing discomfort, and reduced bracing compliance. |
| How long will each pair of AFOs last? | On average, an AFO will last 3-9 months. Some children grow rapidly and will need new AFO's within three months. As the child's growth rate slows, more time can be expected from each AFO. |
| How long will the bar last? | The bars are adjustable and should last approximately more than one pair of AFOs. The bar width should be periodically checked and adjusted to match the child's shoulder width or as recommended by your physician or orthotist. |

8 GLOSSARY


| Term | Meaning |
|--------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Ankle Foot Orthosis (AFO) | An AFO is a brace, usually made of plastic, that is worn on the lower leg and foot to support the ankle, hold the foot and ankle in the correct position |
| Casting | A procedure to keep the clubfoot in place after manipulation. Padding is wrapped around the foot, then plaster or fiberglass is applied, and the cast is molded into place. |
| Clubfoot (talipes equinovarus) | Clubfoot (Talipes equinovarus) is a common birth defect in which the foot is twisted out of shape or position. In clubfoot, the foot appears twisted and can even look as if it's upside down. Despite its appearance, clubfoot itself doesn't cause any discomfort or pain. Treatment is usually successful and includes stretching and casting (Ponseti Method) or stretching and taping (French method). Sometimes, surgery is needed. |
| Dorsiflexion | Dorsiflexion is the flexion of the foot in an upward direction, from the ankle. |
| Equinus | Equinus is a condition in which the upward bending motion of the ankle joint is limited. Someone with equinus lacks the flexibility to bring the top of the foot toward the front of the leg. Equinus can occur in one or both feet. |
| Foot Abduction Brace | Abduction is a motion that pulls a structure or part away from the midline of the body. An abduction brace keeps the foot in proper alignment with the midline. |
| Orthotist | An orthotist is a healthcare professional who makes and fits braces and splints (orthoses) for people who need added support for body parts that have been weakened by injury, disease, or disorders of the nerves, muscles, or bones |
| Ponseti Method | The Ponseti Method is a non-surgical clubfoot treatment process developed by Dr. Ignacio Ponseti. The method utilizes gentle manipulation and casting to gradually move the foot into the correct position. |
| Pressure Saddles | A pressure saddle is a piece of soft, flexible rubber that attaches to the middle strap of the AFO to provide additional relief from the pressure of the strap. |
| Tenotomy | An Achilles tendon release, also known as tenotomy, is a surgical procedure that involves cutting through the Achilles tendon to allow for final dorsiflexion of the foot. The procedure is used to correct the equinus deformity seen in clubfoot. |

9 ICON LEGEND

| Icon | Meaning | Icon | Meaning |
|-----------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|---------------------------------|
|  | Ankle-Foot Orthotic (AFO) |  | Caution |
|  | Authorized Representative in the European Community |  | Date of manufacture |
|  | Manufacturer |  | Importer of record |
|  | Single Patient, Multiple Use |  | Non-sterile |
|  | Indicates Item is a Medical Device |  | Catalog number |
|  | CE Mark (in compliance with Regulation (EU) 2017/745 of the European Parliament and of the council of 5 April 2017 on Medical Devices) |  | Batch code |
|  | Consult Instructions for Use |  | Electronic Instructions for Use |

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MD Orthopaedics, Inc.
604 North Parkway Street
Wayland, IA 52654, САЩ
www.mdorthopaedics.com

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Copyright © 2022 MD Orthopaedics, Inc.