## INSTRUCTION FOR USE

**Mitchell Ponseti® Bracing System**  
**Ponseti® Ankle Foot Orthotic**  
**Ponseti® Plus Abduction Bar**  
**Published: 2023**

Document MD-175  
Version 2  
January 2024


**MD Orthopaedics, Inc.**  
604 North Parkway Street  
Wayland, IA 52654 USA  
[www.mdorthopaedics.com](http://www.mdorthopaedics.com)

MedEnvoy Switzerland  
Gottthardstrasse 28  
6302 Zug  
Switzerland

Emergo Europe B.V.  
Westervoortsedijk 60  
6827 AT Arnhem  
The Netherlands

Emergo Consulting (UK) Limited  
Compass House, Vision Park  
Hilston  
c/o Cr260-UL International  
Cambridge  
England, United Kingdom  
CB24 9BZ

MedEnvoy  
Prinses Margrietplantsoen 33  
Suite 123  
2595 AM The Hague  
The Netherlands

## Table of Contents

English .....	3
Français (Canada) .....	6
Español (EE. UU./Latinoamérica) .....	9
Deutsch .....	12
Español (UE).....	15
Français (UE) .....	18
Italiano.....	21
Nederlands .....	24
Português (EU) .....	27
Ελληνικά.....	30
Eesti.....	33
Magyar.....	36
Lietuviškai.....	39
Latviski.....	42
Română .....	45
Türkçe .....	48
Dansk .....	51
Suomi.....	54
Íslenska.....	57
Norsk.....	60
Svenska.....	63
Български.....	66
Hrvatski .....	69
Čeština .....	72
Polski.....	75
Slovenčina.....	78
Slovenščina .....	81
简体中文.....	84
日本語.....	87
العربية.....	90
עברית.....	93
Afrikaans .....	96
IsiXhosa .....	99
IsiZulu .....	102
Northern Sotho.....	105

## **Intended Use**

The **Ponseti® Plus Abduction Bar** is a Class I, non sterile medical device intended for use only on children whose clubfoot is being corrected with the Ponseti Method of treatment, and only as prescribed by a physician fully trained in the Ponseti Method. This **Ponseti® Plus Abduction Bar** should be used in conjunction with the **Mitchell Ponseti® Ankle Foot Orthosis (AFO)**.

## **Indications**

The **Mitchell Ponseti® Ankle Foot Orthosis (AFO)** with the **Ponseti® Plus Abduction Bar**, is an ankle-foot orthosis system, also known as **Mitchell Ponseti® Bracing System** and is intended for use in the Ponseti Method of treatment for the correction of Congenital Talipes Equinovarus (CTEV) in infants and children less than or equal to eight years of age. The brace is intended to prevent relapse (equinus and varus deformity of the heel).

The **Ponseti® Plus Abduction Bar** is available in 10, or 15 degrees (a 0 degree quick clip is offered as a component) of dorsiflexion and should be used with the **Mitchell Ponseti® Ankle Foot Orthosis**.

The Mitchell Ponseti® Bracing System is to be prescribed by a physician or advanced practice provider trained in the Ponseti Method of clubfoot treatment and is typically worn for 23 hours a day for three months after casting and then during night-time and nap time for 4-5 years or until the physician or advanced practice provider indicate the device is no longer needed.

A clinical orthotist, also trained in the Ponseti Method, may assist the physician or advanced practice provider in determining the correct size of the AFO for the child. The orthotist may also train the child's parents and care givers in the right way to put the AFO onto the child's foot, what cues may indicate poor fit or other issues, and how to determine it is time to move to the next size.

Other members of the clinical team trained in the Ponseti Method – nurses, physician assistants, or other licensed healthcare professionals – may also train parents on the proper use of the Mitchell Ponseti® Bracing System and make any necessary adjustments in fit or sizing. The child's parents, family or other designated caretaker are responsible for adhering to the prescribed Ponseti Method treatment regimen.

## **Clinical Benefits**

The MD Orthopaedics Mitchell Ponseti® Bracing System retains the position of the anatomically corrected clubfoot deformity for up to 4-5 years (or until the physician or advanced practice provider indicate the device is no longer needed) following successful manipulation and casting as described by Dr. Ignacio Ponseti (the Ponseti Method), with a relatively low risk of relapse/recurrence and/or surgical intervention.

## **Manufacturing Information**

Instructions for Use (IFU) manuals are available in English and other languages at [www.mdorthopaedicsifu.com](http://www.mdorthopaedicsifu.com) (Regulation (EU) 207/2012).

MD Orthopaedics, Inc.

604 North Parkway Street

Wayland, IA 52654 USA

[www.mdorthopaedics.com](http://www.mdorthopaedics.com)

1-877-766-7384

The manufacturer is not liable for cases of material damage or personal injury caused by incorrect handling or non-compliance with instructions. In such cases, the warranty will be voided. Warranty of the Mitchell Ponseti® Bracing System is one year under normal use. Normal use is defined as a single user following intended use.

## **Cleaning Instructions**

### **Ponseti® Plus Abduction Bar**

Cleaned by using a mild detergent and damp cloth. Do not submerge.

### **Mitchell Ponseti® Ankle Foot Orthosis (AFO)**


Hand or machine wash using cool water, and mild detergent. A gentle washer cycle may be used. Air dry, do not place in dryer. Do not use hot water when washing.

## **Preparation**

The Mitchell Ponseti® Bracing System is shipped at room temperature and requires no special handling during transport or storage.

## **Instructions**

### **Mitchell Ponseti AFO**

Step 1

Step 2

Step 3


Step 4

Step 5

Step 6

1. Insert the child's foot into the open AFO, with the buckles on the inside (facing the other foot).
2. While holding the heel in place, place the tongue across the front of the leg and top of the foot, as shown.
3. Secure the middle strap over the tongue and buckle it firmly in place. The hole in the tongue should still be visible just above the strap.
4. Fasten the toe strap, making sure that the child's foot is flat and the toes are straight.
5. Fasten the top strap, making sure it is just tight enough to keep the AFO in place.
6. Check the hole at the back of the AFO for proper placement of the heel.

### Ponseti® Plus Abduction Bar

1. It is important to attend all scheduled provider appointments to ensure that shoulder measurements and bar width adjustments are completed properly. For sizing information email [info@mdorthopaedics.com](mailto:info@mdorthopaedics.com).
2. A 5/64" allen (hex) key is provided for bar lock screws. Turn right to tighten until bars are locked into place and left to loosen.
3. The width of the bar should match the shoulder width of the child once locked into place.
4. A #2 Phillips head screw driver fits the center screw of the heel cap.
5. Insert clip into sole of AFOs until you hear a "click". Verify that the Quick Clips are engaged by pulling on the AFO to make sure it does not disconnect from the bar. To disconnect the AFO from the Bar push firmly on "PUSH" while pulling back on the bar until it slides free.

### **Inherent Risks & Contraindications**

**Never use an AFO or bar that is damaged, broken or not functioning properly.**

Modification is not recommended for the Mitchell Ponseti® Bracing System.

Protect your child, yourself, and furniture using a Ponseti® Bar Cover.

### **Product Components and Accessories**

The Mitchell Ponseti® Ankle Foot Orthotic should be used with the following original accessories and components only:

- o Ponseti® Abduction Bar
- o Ponseti® Bar Covers
- o 0° Quick Clip
- o Ponseti® Plus Abduction Bar
- o Ponseti® Pressure Saddles
- o Mitchell Ponseti® Move Bar
- o Ponseti® Socks

MD Orthopaedics does not provide medical treatment, advice or recommendations about the risks and benefits of medical treatment, including treatment that involves the use of MD Orthopaedics products. This information should be provided solely by the physician or other qualified health care provider treating your child. If you have questions about your child's treatment, it is important for you to discuss those questions with the appropriate health care provider.

### **Safety Information Regarding Use**

- o Follow the qualified health care provider and physician instructions or treatment plan.
- o Consult your physician or clinical team immediately if the child is experiencing pain, blisters, sores, or if the AFO is visibly slipping on the child's foot.
- o Socks should be worn when wearing the Mitchell Ponseti® Ankle Foot Orthotic.
- o Children should not stand or walk in the Mitchell Ponseti® Bracing System.

If a user and/or patient experiences any serious incident that has occurred in relation to the device, it should be reported to the MD Orthopaedics, Inc. and the competent authority of the EU and EEA Member State in which the user and/or patient is established.

### Safe Disposal

The Mitchell Ponseti® Bracing System can be disposed of in the regular trash. None of the components of the device are recyclable.

### Reusing the Device

All parts were designed to be used by a single patient under normal use conditions and as prescribed by a physician. Normal use includes the daily application of the Mitchell Ponseti® Bracing System and use up to 23 hours per day, or as prescribed.

### Symbol Legend

Symbol	Meaning	Symbol	Meaning
	Ankle-Foot Orthotic (AFO)		Caution
	Authorized Representative in the European Community		Date of manufacture
	Manufacturer		Importer of record
	Single Patient, Multiple Use		Non-sterile
	Indicates Item is a Medical Device		Catalog number
	CE Mark (in compliance with Regulation (EU) 2017/745 of the European Parliament and of the council of 5 April 2017 on Medical Devices)		Batch code
	Information considered important, but not hazard-related		Electronic Instructions for Use

**MD Orthopaedics, Inc.**  
604 North Parkway Street  
Wayland, IA 52654 USA  
www.mdorthopaedics.com

MedEnvoy Switzerland  
Gotthardstrasse 28  
6302 Zug  
Switzerland

Emergo Europe B.V.  
Westervoortsedijk 60  
6827 AT Arnhem  
The Netherlands

Emergo Consulting (UK) Limited  
Compass House, Vision Park  
Hilston  
c/o Cr3600-UL International  
Cambridge  
England, United Kingdom  
CB24 9BZ

MedEnvoy  
Prinses Margrietplantsoen 33  
Suite 123  
2595 AM The Hague  
The Netherlands

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