**MANUFACTURER:**  
ORTHOPEDIATRICS  
2850 Frontier Drive  
Warsaw, IN 46582 USA  
Tel: 574-268-6379

**AUSTRALIAN SPONSOR:**  
Orthopediatrics AUS Pty Ltd  
15 Orchard Road, Suite 1  
Brookvale NSW 2100

**EUROPEAN AUTHORIZED  
REPRESENTATIVE**  
Emergo Europe  
Westervoortsedijk 60  
6827 AT Arnhem The Netherlands

CE Mark, UKCA Mark, and Single Use symbols are valid only if also printed on the product label or device itself

[www.orthopediatrics.com/symbol\\_glossary](http://www.orthopediatrics.com/symbol_glossary)

[www.orthopediatrics.com/ifu](http://www.orthopediatrics.com/ifu)

[IMPORTANT MEDICAL INFORMATION – PACKAGING INSERT](#)

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The Dynamic Femur Fracture (DF2) brace is an individually fitted brace comprised of several components including (all or some of) the following: pelvic girdle, hip hinge, knee hinge, outer sleeve, inner sleeve, tongues, and foot plate.

**MATERIALS**

The DF2 brace is made from various types of metal and plastic materials. The materials include low density polyethylene, polypropylene, stainless steel, aluminum, nickel plated brass, nylon, pressure sensitive adhesive, polyethylene terephthalate, crosslinked polyolefin foam, polyvinyl chloride, hook-and-loop fasteners, foam laminate, polyester, and nickel.

**INDICATIONS FOR USE**

The DF2 brace is intended for femur fracture fixation and post-operative stabilization in pediatric patients from approximately 6 months to 5 years of age by providing immobilization of the femur, knee, and hip.

**CONTRA-INDICATIONS**

Do not use the DF2 brace in the presence of any contraindication. Contraindications include but are not limited to:

- Material sensitivity documented or suspected.
- Inability to follow treatment protocol and care instructions.

**WARNINGS**

- Federal Law restricts the DF2 brace to sale by, or on the order of, a qualified physician or healthcare provider and is only to be used as prescribed.
- Do not use if the fracture cannot be fully encapsulated in the brace.
- Do not remove subcomponents of the DF2 brace as they may present a choking hazard. Ensure proper assembly to prevent loose components around small children.
- Do not use the DF2 brace without a stockinette underneath the brace. This will prevent skin breakdown and improve overall hygiene during patient's time in the brace.
- It is necessary that the brace be worn for 24 hours per day or as directed by your physician or healthcare provider.
- Incorrect use of the DF2 brace may result in poor fracture healing, malunion, or nonunion.
- Extra care should be used with brace when used with any open wounds and/or infections.
- The DF2 brace a single-patient reusable device and shall not be used on multiple patients.

**POSSIBLE ADVERSE EFFECTS**

Consult your physician or healthcare provider immediately if the patient is experiencing any of the following:

- Skin irritation
- Pressure sores
- Pain, discomfort, or abnormal sensations due to the presence of the device
- Increased swelling
- Signs of infection
- Brace failure

These risks should be explained to the patient. Patient should contact a health professional in relation to the operation of the device.

**MRI SAFETY INFORMATION**

The DF2 brace is MR unsafe and should remain outside the MRI scanner room.

## CLEANING INSTRUCTIONS

1. Do not fully remove the DF2 Brace for cleaning unless directed by a physician or healthcare provider. Hook-and-loop straps may be undone to access the inside of the brace for cleaning purposes.
2. Stockinette changes will be directed by your physician or healthcare provider. If the stockinette becomes soiled, reach out to your healthcare provider for changing instructions.
3. The DF2 brace may be cleaned with a washcloth and warm water.
4. The DF2 brace may be disinfected with rubbing alcohol or mild antibacterial wipes.
5. Harsh cleaners, such as bleach, are not recommended as the brace is in close contact with the skin.
6. Do not put the DF2 brace into a mechanical dryer. Air dry only.
7. Do not submerge the DF2 brace in water or cleaning liquids.

## REUSE LIFE

All parts were designed to be used by a single patient under normal use conditions and as prescribed by a physician or healthcare professional. Normal use includes full time wearing of the DF2 brace as prescribed.

## STERILIZATION INFORMATION

The DF2 brace cannot be sterilized via steam, gamma, ETO, or x-ray. The brace does not require sterilization before use.

## STORAGE AND HANDLING

- The DF2 brace requires no special handling during transport or storage.

## IMPORTANT STATEMENT

- In case of any serious incident that has occurred in relation to the device, the user and/or patient should report the incident to OrthoPediatrics and the EU competent authority of the Member State or Therapeutic Goods Administration (<https://www.tga.gov.au/>) in which the user and/or patient is established.
- Please see the User Manual for additional information. An electronic copy can be found at <https://www.orthopediatrics.com/>. To obtain a copy of the User Manual please call the Customer Service Group at OrthoPediatrics at 574-268-6379.
- DF2 braces should be checked for defects before use and to verify the appropriate sizing.
- The Dynamic Femur Fracture (DF2) Brace and accessories can be disposed of in accordance with applicable local, state, and federal regulations. None of the components of the device are recyclable.
- Only a qualified physician or healthcare professional may modify the DF2 brace.

July 2024

OrthoPediatrics is a trademark, registered in the United States, of OrthoPediatrics Corp.

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OrthoPediatrics je blagovna znamka, registrirana v Združenih državah Amerike, OrthoPediatrics Corp.