**Mitchell Ponseti®**

MD Orthopaedics

## Instructions for Use

**MITCHELL PONSETI® Ankle Foot Orthosis (AFO)**

**PONSETI® Abduction Bar**

**MITCHELL PONSETI® Move Bar**

**MD Orthopaedics, Inc.**  
604 North Parkway Street  
Wayland, IA 52654 USA  
[www.mdorthopaedics.com](http://www.mdorthopaedics.com)

Document MD-124.3  
Version 10  
March 2024

Emergo Europe B.V.  
Westervoortsedijk 60  
6827 AT Arnhem  
The Netherlands

MedEnvoy  
Prinses Margrietplantsoen 33  
Suite 123  
2595 AM The Hague  
The Netherlands

MedEnvoy Switzerland  
Gotthardstrasse 28  
6302 Zug  
Switzerland

## Table of Contents

DISCLAIMER .....	1
TRADEMARKS .....	1
<b>1 PREFACE .....</b>	<b>1</b>
1.1 Intended Purpose .....	1
1.2 Target Patient Group and Intended Users .....	2
1.3 Clinical Benefits .....	2
1.4 Explanation of Safety Warnings .....	2
1.5 Retention Instructions & Warranty.....	2
1.5.1 Technical Life Span .....	2
1.5.2 Returns/Refunds .....	3
1.6 Obtaining Documentation and Information .....	3
1.6.1 Ordering Documentation .....	3
1.6.2 Other languages .....	3
1.6.3 Documentation Feedback .....	3
1.6.4 Support and service .....	3
<b>2 DESCRIPTION OF THE PRODUCT .....</b>	<b>3</b>
<b>3 Clubfoot Treatment Overview .....</b>	<b>5</b>
3.1 MITCHELL PONSETI® Ankle Foot Orthotic .....	6
3.2 MITCHELL PONSETI® AFO .....	7
3.3 PONSETI® Abduction Bar .....	8
3.4 Mitchell Ponseti® Move Bar .....	11
<b>4 SAFETY INSTRUCTIONS .....</b>	<b>13</b>
4.1 Safe Disposal .....	14
4.2 Potential Health Consequences .....	14
<b>5 STORAGE AND TRANSPORT .....</b>	<b>14</b>
5.1 How to Transport and Store the Product .....	14
<b>6 MAINTENANCE .....</b>	<b>14</b>
6.1 Reusing the Device .....	14
6.1.1 Cleaning the device .....	14
6.2 How to Inspect the Product .....	14
<b>7 TROUBLESHOOTING .....</b>	<b>15</b>
7.1 How to Identify and Solve Problems .....	15
7.2 Frequently Asked Questions .....	15
<b>8 GLOSSARY .....</b>	<b>16</b>
<b>9 SYMBOL LEGEND .....</b>	<b>17</b>

## DISCLAIMER

MD Orthopaedics makes no representations or warranties with respect to this manual and, to the maximum extent permitted by law, expressly limits its liability for breach of any warranty that may be implied to the replacement of this manual with another. Furthermore, MD Orthopaedics reserves the right to revise this publication at any time without incurring an obligation to notify any person of the revision.

As the designer and manufacturer of products, MD Orthopaedics does not provide medical treatment and/or advice. Information about the risks and benefits of medical treatment, including treatment that involves the use of MD Orthopaedics products and available alternatives is provided solely by physicians or other healthcare providers. If you have questions about your child's treatment, it is important to discuss those questions with the appropriate healthcare provider.

MD Orthopaedics does not, and cannot, warrant certain results from the use of its products.

The information provided in this documentation contains general descriptions and/or technical characteristics of the performance of the products contained herein. This documentation is not intended as a substitute for and is not to be used for determining suitability or reliability of these products for specific user applications. It is the duty of any such user to perform the appropriate and complete risk analysis, evaluation, and testing of the products with respect to the relevant specific application or use thereof. Neither MD Orthopaedics nor any of its affiliates or subsidiaries shall be responsible or liable for misuse of the information that is contained herein. If you have any suggestions for improvements or amendments or have found errors in this publication, please notify us.

All pertinent state, regional, and local safety regulations must be observed when using this product. For reasons of safety and to help ensure compliance with documented system data, only the manufacturer shall perform repairs to components.

When devices are used for applications with technical safety requirements, the relevant instructions must be followed. Failure to observe this information can result in injury or equipment damage.

MD Orthopaedics has made every effort to prevent that the labelling, instructions for use, making available, putting into service and advertising of devices uses text, names, trademarks, pictures and figurative or other signs that may mislead the user or the patient regarding the device's intended purpose, safety, and performance.

Copyright © 2023 by MD Orthopaedics, Inc.

All rights reserved. No part of this publication may be reproduced, distributed, or transmitted in any form or by any means, including photocopying, recording, or other electronic or mechanical methods, without the prior written permission of the publisher. For permission requests, write to the publisher at the address below.

MD Orthopaedics, Inc.  
604 North Parkway Street  
Wayland, IA 52654  
USA  
1-877-766-7384  
[www.mdorthopaedics.com](http://www.mdorthopaedics.com)

## TRADEMARKS

**Mitchell Ponseti®** and **Ponseti®** are registered trademarks of MD Orthopaedics.

MD Orthopaedics has made every effort to supply trademark information about company names, products and services mentioned in this manual. Trademarks shown below were derived from various sources. All trademarks are the property of their respective owners.

General Notice: Some product names used in this manual are used for identification purposes only and may be trademarks of their respective companies.

## 1 PREFACE

### 1.1 Intended Purpose

The **MITCHELL PONSETI® Ankle Foot Orthosis (AFO)** with the **PONSETI® Abduction Bar** or the **MITCHELL PONSETI® Move Bar** is known as **MITCHELL PONSETI® Brace**, intended for use in the Ponseti method of treatment for the correction of Congenital Talipes Equinovarus (CTEV) in infants and children less than or equal to eight years of age. The brace is intended to prevent relapse (equinus and varus deformity of the heel).

The **MITCHELL PONSETI® Brace** is to be prescribed by a physician or healthcare provider trained in the Ponseti Method of clubfoot treatment and is typically worn for 23 hours a day for three months after casting and then during night-time and nap time for 4-5 years or until the physician or healthcare provider indicate the device is no longer needed.

A clinical orthotist, also trained in the Ponseti method, may assist the physician or healthcare provider in determining the correct size of the AFO for the child. The orthotist may also train the child's parents and families in the right way to put the AFO onto the child's foot, what cues may indicate poor fit or other issues, and how to determine it is time to move to the next size.

Other members of the clinical team trained in the Ponseti method – nurses, physician assistants, or other licensed healthcare professionals – may also train parents on the proper use of the **MITCHELL PONSETI® Brace** and make any necessary adjustments in fit or sizing.

The child's parents or family or other designated caretaker are responsible for adhering to the prescribed Ponseti method treatment regimen of 23 hours a day for three months, then during sleep and during naptime for the next 4-5 years or until the physician indicates the device is no longer needed.

## 1.2 Target Patient Group and Intended Users

The MD Orthopaedics **MITCHELL PONSETI® Ankle Foot Orthosis (AFO)** with the **PONSETI® Abduction Bar** or **MITCHELL PONSETI® Move Bar** is designed to be worn by a child ( $\leq 8$ ) years of age with clubfoot after the casting portion of treatment is complete.

This document is intended for the physician, orthotist, clinical team, parents, and families of children receiving treatment with the **MITCHELL PONSETI® Ankle Foot Orthosis (AFO)** with the **PONSETI® Abduction Bar** or **MITCHELL PONSETI® Move Bar**.

## 1.3 Clinical Benefits

The MD Orthopaedics **MITCHELL PONSETI® Brace** (that includes the **MITCHELL PONSETI® Ankle Foot Orthosis (AFO)** and either the **PONSETI® Abduction Bar** or the **MITCHELL PONSETI® Move Bar**) retains the position of the anatomically corrected clubfoot deformity for up to 4-5 years following successful manipulation and casting as described by Dr. Ignacio Ponseti (the Ponseti Method), with a relatively low risk of relapse/recurrence and/or surgical intervention.

## 1.4 Explanation of Safety Warnings

 Caution indicates a hazard with a low level of risk which, if not avoided, could result in minor or moderate injury.

 Indicates information considered important, but not hazard-related.

## 1.5 Retention Instructions & Warranty

Read and understand this manual and its safety instructions before using this product. Failure to do so can result in injury.

The product shall only be used by persons who have fully read and understand the contents of this user manual.

Ensure that each person who uses the product has read these warnings and instructions and follows them.

The manufacturer is not liable for cases of material damage or personal injury caused by incorrect handling or non-compliance with the safety instructions. In such cases, the warranty will be voided. Warranty of the **Mitchell Ponseti® Brace** is one year under normal use, and is single patient multiple use. Normal use is defined by Normal Daily Use for 365 days/year.

This policy/warranty is non-transferable and covers normal wear and tear only. Altering the product and after-market modifications will void the warranty. Please contact us if you experience problems not outlined in this policy. Atypical matters will be handled on a case-by-case basis.

### 1.5.1 Technical Life Span

On average, each AFO will last 3-9 months, depending on the child's growth rate. As the child's growth rate slows, more time can be expected.

The bars are adjustable and may last more than one pair of AFOs. The bar width should be periodically checked and adjusted to match the child's shoulder width or as recommended by your physician, orthotist or healthcare provider.

### 1.5.2 Returns/Refunds

If the merchandise you received is defective or not as you ordered, please contact customer service or your point of purchase for quick resolution.

For information on returning product and refunds, please contact

MD Orthopaedics, Inc.

604 North Parkway Street

Wayland, IA 52654

USA

[www.mdorthopaedics.com](http://www.mdorthopaedics.com)

1-877-766-7384 or e-mail

[info@mdorthopaedics.com](mailto:info@mdorthopaedics.com)

No refunds or exchanges will be issued for Custom Color AFOs; all sales are final.

## 1.6 Obtaining Documentation and Information

### 1.6.1 Ordering Documentation

Additional documentation, user instructions and technical information can be ordered by calling MD Orthopaedics at 1-877-766-7384 or e-mailing [info@mdorthopaedics.com](mailto:info@mdorthopaedics.com).

### 1.6.2 Other languages

Instructions for Use (IFU) manuals are available in other languages at [www.mdorthopaedicsifu.com](http://www.mdorthopaedicsifu.com) (Regulation (EU) 207/2012).

### 1.6.3 Documentation Feedback

If you are reading MD Orthopaedics' product documentation on the internet, any comments can be submitted at [www.mdorthopaedics.com](http://www.mdorthopaedics.com) (Regulation (EU) 207/2012) or comments can also be sent to [info@mdorthopaedics.com](mailto:info@mdorthopaedics.com). We appreciate your comments.

### 1.6.4 Support and service

For questions, information, technical assistance or to order user instructions, please contact:

MD Orthopaedics, Inc.

604 North Parkway Street

Wayland, IA 52654

USA

[www.mdorthopaedics.com](http://www.mdorthopaedics.com)

1-877-766-7384 or e-mail [info@mdorthopaedics.com](mailto:info@mdorthopaedics.com)

## 2 DESCRIPTION OF THE PRODUCT

The **MITCHELL PONSETI® Brace** is a Class I medical device intended for use only on children whose clubfoot is being corrected with the Ponseti Method of treatment, and only as prescribed by a physician or healthcare provider fully trained in the Ponseti Method.

**Preemie Bar:** Size P6-0 or P5-0 AFOs are permanently attached to a Preemie Bar, which is pre-set to 10° dorsiflexion and 60° abduction. The bar adjusts in width from 15.0-20.0 cm, in 1-cm increments. The Preemie Bar can be cut down as small as 10 cm, upon request.

The **MITCHELL PONSETI® Ankle Foot Orthosis** shall be used with the following original accessories and components only:

- **PONSETI® Adjustable Bar:** Adjustable bar that easily clips on the AFO with the quick clip system.
  - Available with either 10 or 15 degrees of dorsiflexion
  - Available in Extra Short, Short, and Long
- **MITCHELL PONSETI® Move Bar:** A dynamic bar that allows flexion and extension of each leg while maintaining external rotation of the feet.
  - Available with either 10 or 15 degrees of dorsiflexion
  - Available in Short, Medium, and Long
- **Pressure Saddle:** This soft, comfortable pressure saddle attaches to the middle strap of the **MITCHELL PONSETI® AFO** to provide additional relief of pressure areas.
  - Available in size small and large. Size small is recommended for use with AFO Size 0000-1; Size large is recommended for use with Size 2-12.
- **MITCHELL PONSETI® Socks** are a single layer sock designed to be used with AFOs using a soft, natural antibacterial, moisture wicking bamboo material that have silicone grips and a high-visibility heel.
  - Available in sizes 0 through 12.
- **Bar Cover:** This soft bar cover provides protection and comfort for your child, and is available in size small (5") and large (8"), in pink, light blue, or dark blue.

The **MITCHELL PONSETI® Ankle Foot Orthosis** is available in four styles:

The **MITCHELL PONSETI® Standard** AFO uses the Ponseti Method of clubfoot correction. Our patented system features a soft lining for comfort and compliance, soft synthetic leather straps. The location of the heel can easily be seen through two viewing holes in the back of the AFO. The AFO is worn attached to the **PONSETI® Abduction Bar** or **MITCHELL PONSETI® Move Bar** in accordance with the standard Ponseti bracing protocol. Available in sizes preemie through 12.

The **MITCHELL PONSETI® Plantar Flexion Stop** AFO limits plantar flexion with rigid dorsiflexion to maintain the prescribed position of the foot. Recommended for the complex clubfoot. This device is also good for children with hyper flexible feet as it ensures the foot is always held in a neutral or dorsiflexed position. The AFO is worn attached to the **PONSETI® Abduction Bar** or **MITCHELL PONSETI® Move Bar** in accordance with the standard Ponseti bracing protocol. Available in sizes 00 through 12. The **MITCHELL PONSETI® Move Bar** is compatible with sizes 00 through 5.

The **MITCHELL PONSETI® Toe Stilt** AFO was specifically designed by recommendation of Dr. Ponseti to maintain clinically desirable dorsiflexion by stretching the Achilles tendon and helping to keep the foot flexible. This is achieved by the placement of a built in wedge attached to the front of the AFO which keeps both feet in a dorsiflexed position of 10 degrees when standing. The AFO is worn attached to the **PONSETI® Abduction Bar** or **MITCHELL PONSETI® Move Bar** in accordance with the standard Ponseti bracing protocol. The patient must be able to stand independently for the Toe Stilt AFO to have any effect when used without the **PONSETI® Abduction Bar** or **MITCHELL PONSETI® Move Bar**. Available in sizes 2 through 11.


The **MITCHELL PONSETI® Plantar Flexion Stop/Toe Stilt** is a dual purpose AFO that supports the Ponseti Method of treating clubfoot. The AFO is worn attached to the **PONSETI® Abduction Bar** or **MITCHELL PONSETI® Move Bar** in accordance with the standard Ponsetibracing protocol. The Plantar Flexion Stop (PFS) is designed to structurally limit any plantar flexion and work particularly well in clubfoot cases with hypermobile or corrected atypical/complex deformity. The PFS is mounted on the standard AFO during the assembly process but is designed to be custom fit as necessary by the orthotist to meet the dynamic needs of each clubfoot child. The Toe Stilt encourages lower calf, Achilles tendon and plantar fascia stretch in the toddler who can stand independently without the **PONSETI® Abduction Bar** or **MITCHELL PONSETI® Move Bar**. The Toe Stilt holds both feet at 10 degrees dorsiflexion when standing. Available in sizes 2 through 11. The **MITCHELL PONSETI® Move Bar** is compatible with sizes 00 through 5.

**!** Education on proper brace application and troubleshooting for skin issues is critical to ensure your child's clubfoot correction is effective. DO NOT USE this device without first being trained by your physician or clinical team. Contact Customer Support for assistance in locating training resources.

**!** All products may not be available in your region, please check with your local distributor for more information.

### 3 Clubfoot Treatment Overview

# Clubfoot Journey™

This info graphic provides that such medical information is merely information - not advice. If users need medical advice, they should consult a doctor or other appropriate medical professional. No warranties are given in relationship to the medical information supplied and no liability will be accrued to MD Orthopaedics in the event the user suffers loss as a result of reliance upon the information.

*This infographic provides information — not advice. If you need medical advice, please consult a doctor or other appropriate medical professional.*

*No warranties are given in relationship to the medical information supplied and no liability will be accrued to MD Orthopaedics in the event the user suffers loss because of reliance upon the information provided in this infographic.*

#### Technical Characteristics


The **MITCHELL PONSETI® AFO** was designed with a soft, contoured elastomeric liner to cushion and hold the child's foot in place. Straps and a body made from a ductile synthetic suede that conforms to the shape of the foot for added comfort and stability. The AFO has the openings on the side to add ventilation and keep the foot cool, and viewing holes to help gauge proper placement of the heel in the AFO.

 Socks are **required** to eliminate skin contact with product.

The Quick Clips on the **PONSETI® Abduction Bar** and the **MITCHELL PONSETI® Move Bar** have been factory set to 60 degrees, the most common bilateral clubfoot setting. This configuration may not correspond to your physician's preferred clubfoot setting. Please discuss with your physician, orthotist, clinical team or healthcare provider before using this device.

### 3.1 MITCHELL PONSETI® Ankle Foot Orthotic

- a. Liner
- b. Tongue
- c. Straps
- d. Buckles
- e. Sole


MD Orthopaedics' AFO Sizing Chart

MITCHELL PONSETI® AFO Size	(A) Foot Length (cm)	(B) Foot Width (cm)
P6-0	5.5 to 6.5	2.5 to 4.0
P5-0	5.8 to 6.9	3.0 to 4.5
0000	6.0 to 7.3	3.3 to 4.8
000	6.3 to 7.8	3.0 to 5.5
00	6.8 to 8.3	3.0 to 5.5
0	7.5 to 9.0	3.0 to 6.0
1	8.0 to 9.5	3.5 to 6.0
2	9.0 to 10.5	3.8 to 6.5
3	10.0 to 11.5	4.5 to 7.0
4	11.0 to 12.5	4.7 to 7.0
5	12.0 to 13.5	4.8 to 7.0
6	13.0 to 14.5	5.0 to 7.5
7	14.0 to 15.5	5.5 to 8.0
8	15.0 to 16.5	6.7 to 8.2
9	16.0 to 17.5	7.0 to 8.5
10	17.0 to 18.5	7.4 to 8.9
11	18.0 to 19.5	7.7 to 9.2
12	19.0 to 20.5	8.0 to 9.5

Legacy European AFO Sizing Chart

MITCHELL PONSETI® AFO Size	(A) Foot Length (mm)
P6-0 (Premie 1)	60 to 62
P5-0 (Premie 2)	63 to 66
0000	67 to 70
000	71 to 75
00	76 to 80
0	81 to 87
1	88 to 92
2	93 to 102
3	103 to 112
4	113 to 122
5	123 to 132
6	133 to 142
7	143 to 152
8	153 to 162
9	163 to 172
10	173 to 182
11	183 to 192
12	193 to 202

Length and width measurement locations

**When fitting the AFO it is recommended to size within the sizing parameters of the sizing chart. Please email us for sizing assistance at [info@mdorthopaedics.com](mailto:info@mdorthopaedics.com).**


### 3.2 MITCHELL PONSETI® AFO

Carefully inspect the AFO before each use. Never use an AFO or bar that is damaged, broken or not functioning properly, as these may injure your child and/or render treatment ineffective. Always use socks with the AFO to prevent direct contact with the skin.

 Ensure the area is well lit and free from distraction to ensure you apply the AFO properly. Improper application of the AFO or bar can lead to skin irritation, sores, and blisters.

1. Open the AFO so all the suede straps are clear for inserting the foot. Buckles are designed to be placed to the inside of the foot. They may be positioned to the outside based on physician preference or individual convenience of application.

2. Holding the lower leg, gently slide the foot into place until the heel is secure against the back and bottom of the AFO. Socks that cover the foot and lower leg are required to prevent skin contact with the AFO.


 The exact position of the heel in the AFO will vary depending on the type and severity of the clubfoot. Forcing the heel into position can lead to sores, so check with your physician to verify the proper position.

3. Pull tongue horizontally against the ankle and hold in place with your thumb. Be sure the hole in the tongue is set in the center of the ankle, above the middle strap.

4. Buckle the middle strap securely over the horizontal tongue. Securely buckle the ankle strap but be careful not to overtighten, as that can lead to skin irritation and sores.

5. Look through the heel viewing holes at the back of the AFO to ensure the heel is ideally down and to the back of the AFO. The heel may not initially be able to contact the bottom of the AFO after the tenotomy and last casting but will over time. If you do not see the heel in the viewing holes, readjust the tongue and middle strap so that the heel is properly positioned.

 In some cases, the heel will not touch the bottom of the AFO, and in others the heel will be barely visible through the heel viewing holes. In these challenging cases, the heel will drop over time as long as ankle motion is adequate. Forcing the heel into the AFO and overtightening the middle strap can lead to skin irritation, sores, and reduced bracing compliance.


Once upper straps are tight and heel is securely in place, buckle the toe strap and retighten the other straps if needed.

### 3.3 PONSETI® Abduction Bar

The **PONSETI® Abduction Bar** was designed to be easily adjustable to the shoulder width of the child. The Quick Clips make it easy for the bar to be attached and removed from the AFOs as needed.

- A Quick Clip
- B Adjustable bar
- C Ponseti® Bar Lock
- D Heel Cap
- E Hex Key


PONSETI® Abduction Bar Size	Minimum Width (cm)	Maximum Width (cm)
Extra Short	17.3	25.0
Short	20.0	30.0
Long	23.5	37.8


 The **PONSETI® Abduction Bar**, size Extra Short Bar can be cut down as small as 11 cm, upon request.

The bar should be set so that the width of the bar is equal to the shoulder width of the child. Adjust the length of the bar, measuring from the center screw of the left heel cap to the center screw of the right heel cap, so that it is equal to the shoulder width measurement.


- I** It is important to measure the shoulder width regularly and adjustments should be approved by orthotist. When sizing the Mitchell Ponseti® Abduction Bar, we recommend measuring from the right shoulder girdle to the left as you face the child (as recommended by the Ignacio Ponseti Foundation). See your Ponseti trained provider or your orthotist to make sure your width is correct. Please email us for sizing assistance at [info@mdorthopaedics.com](mailto:info@mdorthopaedics.com).


1. Use the enclosed hex key to loosen set screws in bar lock. Adjust the width of the bar to the child's shoulder width, using the center screw in the heel caps as a guide. Re-tighten the set screws in the bar lock.


2. Using a screwdriver, loosen center screw of heel cap; remove outer screw. Rotate clip to desired degree indicator to set abduction (external rotation) as determined by your physician. Re-insert outer screw into appropriate hole. Tighten both screws.

3. Insert clip into sole of AFOs until you hear a “click”. Verify that the bar lock and Quick Clips are engaged when assembled by pulling on the AFO to make sure it does not disconnect from the bar. To disconnect the AFO from the Bar push firmly on “PUSH” while pulling back on the bar until it slides free. If damage is found discontinue use and contact Customer Service.

### 3.4 Mitchell Ponseti® Move Bar

- A Quick Clip
- B Spring Capsule
- C Mitchell Ponseti® Bar Lock
- D Adjustable Bar
- E Arm
- F Hex Key


MITCHELL PONSETI® Move Bar	Minimum Width (cm)	Maximum Width (cm)
Short	19.2	22.5
Medium	22.4	28
Long	27	37

The **MITCHELL PONSETI® Move Bar** has a dynamic range of 90 degrees. The mild spring tension will be easily overcome by your child's natural positioning. When in 0 degree/neutral position, the feet will be in a 10 or 15 degree dorsiflexion based on healthcare provider's preference.



The bar should be set so that the width of the bar is equal to the measurements of the child's shoulders.

- !** It is important to measure the shoulder width regularly and adjustments should be approved by orthotist. When sizing the **MITCHELL PONSETI® Move Bar**, we recommend measuring from the right shoulder girdle to the left as you face the child (as recommended by the Ignacio Ponseti Foundation). See your Ponseti trained provider or your orthotist to make sure your width is correct. Please email us for sizing assistance at [info@mdorthopaedics.com](mailto:info@mdorthopaedics.com).

1. Use the enclosed hex key to loosen set screw in bar lock. Adjust the width of the bar to the child's shoulder width, using the marking on the bar. Re-tighten the set screw in the bar lock.


2. Using a screwdriver, loosen screw until able to turn the quick clip. Rotate clip to desired degree indicator to set abduction (external rotation) as determined by your physician. Tighten screw. Repeat process with other quick clip.

3. Insert clip into sole of AFOs until you hear an audible “click”. Verify that Quick Clips are engaged when assembled by pulling on the AFO to make sure it does not disconnect from the bar. To disconnect the AFO from the Bar push firmly on “PUSH” while pulling back on the bar until it slides free. If damage is found, discontinue use and contact Customer Service.

**!** When inserting clip into the AFO ensure grooves are engaged.

#### 4 SAFETY INSTRUCTIONS

Use the **MITCHELL PONSETI® Ankle Foot Orthosis (AFO)**, **PONSETI® Abduction Bar** and the **MITCHELL PONSETI® Move Bar** only in accordance with these instructions and as directed by your physician, healthcare provider or other member of the clinical team.

##### **!** Precautions/Cautions and Warnings:

- Read and understand this manual and its safety instructions before using this product. Failure to do so can result in injury.
- Use the **MITCHELL PONSETI® Brace** only in accordance with these instructions and as directed by your physician or other member of the clinical team.
- Consult your physician, clinical team, or healthcare provider immediately if the child is experiencing pain, blisters, sores, visibly slipping or if they are not in good condition or do not correctly fit the child.
- Do not bend the bar or allow children to walk with the bar attached. Bending of the bar may compromise the therapeutic value, lead to breakage and could negatively affect the correction. Walking with the bar attached puts the child at a high fall risk and may cause breakage to the brace which could negatively affect correction.
- Do not over-tighten the adjusting screws.

- Do not use the **MITCHELL PONSETI® Brace** if it is not in good condition or does not correctly fit the child.
- Protect your child, yourself, and your furniture from being hit by the bar when the child is wearing it. It is recommended to pad the bar with a **PONSETI® Bar Cover**.
- The **MITCHELL PONSETI® Brace** shall not be used without prescription by a physician or healthcare provider trained in the Ponseti Method and is to be used only as prescribed.
- The **MITCHELL PONSETI® Brace** shall not be used to treat any other orthopedic condition other than clubfoot unless otherwise prescribed by the physician or healthcare provider.
- Altering the product and/or making after-market modifications will void the warranty and may result in injury.
- Ensure the dressing area is well lit and free from distraction to ensure you apply the AFO properly.
- Improper application of the AFO or bar can lead to skin irritation, sores, and blisters.
- The **MITCHELL PONSETI® Brace** contains small pieces that may present a choking hazard. Ensure proper assembly to prevent loose components around small children.

#### 4.1 Safe Disposal

- The **MITCHELL PONSETI® Brace** and accessories can be disposed of in accordance with applicable local, state and federal regulations. None of the components of the device are recyclable.

#### 4.2 Potential Health Consequences

It is necessary that the brace be worn 23 hours a day for three months and then 12-14 hours a day (naps and nighttime) for 4-5 years or as directed by your physician or healthcare provider. Incorrect use of the AFO may result in relapse of the clubfoot defect or injury to the child's foot. If a user and/or patient experiences any serious incident that has occurred in relation to the device, it should be reported to the MD Orthopaedics, Inc. and the competent authority of the EU and EEA Member State in which the user and/or patient is established.

Please contact MD Orthopaedics if you experience problems not outlined in this policy.

MD Orthopaedics, Inc.

604 North Parkway Street

Wayland, IA 52654 USA

www.mdorthopaedics.com

+1-877-766-7384 or e-mail info@mdorthopaedics.com

## 5 STORAGE AND TRANSPORT

### 5.1 How to Transport and Store the Product

The **MITCHELL PONSETI® Ankle Foot Orthosis**, **PONSETI® Abduction Bar** and the **MITCHELL PONSETI® Move Bar** ship at room temperature and require no special handling during transport or storage.

## 6 MAINTENANCE

### 6.1 Reusing the Device

All parts were designed to be used by a single patient under normal use conditions and as prescribed by a physician.

Normal use includes the daily application of the AFO and **MITCHELL PONSETI® Brace** and use up to 23 hours per day, or as prescribed.

#### 6.1.1 Cleaning the device

To clean the device:

1. Hand or machine-wash the AFO using cool water, a gentle cycle and mild detergent.
2. Let the AFO air dry completely before use.
3. It is not necessary to wash before using.
4. **PONSETI® Abduction Bar** and **MITCHELL PONSETI® Move Bar** can be properly cleaned by using a mild detergent and damp cloth.

 Do not use hot water when washing the **MITCHELL PONSETI® Ankle Foot Orthosis**. This may cause the synthetic suede straps to warp, shrink, separate or delaminate. Do not put the AFO into a mechanical dryer. Air dry only.

 **PONSETI® Abduction Bar** and **MITCHELL PONSETI® Move Bar** should not be submerged.

### 6.2 How to Inspect the Product

The **MITCHELL PONSETI® Ankle Foot Orthosis (AFO)** and **PONSETI® Abduction Bar** should be inspected every time for proper fit, manufacturing defect, and damage prior to the device being placed onto the child's feet.

 Carefully inspect the **MITCHELL PONSETI® Ankle Foot Orthosis (AFO)**, **PONSETI® Abduction Bar** and the **MITCHELL PONSETI® Move Bar** before each use. Never use an AFO or bar that is damaged, broken or not functioning properly, as this may injure your child and/or render treatment ineffective. Contact Customer Service for Support, if needed.

## 7 TROUBLESHOOTING

### 7.1 How to Identify and Solve Problems

Hazard	Cause	Solution
Child develops blisters or other damage to the skin	Improper fit and/or use	Contact your physician or healthcare provider
Broken AFO	Manufacturing defect or damaged via normal use	Contact your point of purchase
Broken AFO	Damaged via normal use	Contact your point of purchase
Broken Bar, Quick Clip, Heel Cap or Bar Lock	Damaged via normal use	Contact your point of purchase
Incorrect Size (within 7 days of receipt)	Sizing Error	Contact your point of purchase

### 7.2 Frequently Asked Questions

Question	Answer
Do the buckles go to the inside of the foot or to the outside of the foot?	The AFO was designed with the buckles set to the inside of the foot. This design feature is to aid parents in the donning and doffing of the AFOs. Ponseti experts do allow buckles to the outside of the foot based on professional preference or customer convenience. It will not alter the performance or risk of injury by orienting the buckles outward.
I just received my order, and the AFOs seem too big/too small. Is there a different size that would work better?	Contact your point of purchase
The child's heel is not down in the AFO. Is something wrong?	If the heel is not visible in the viewing holes, check with your physician to ensure the clubfoot is fully corrected. If the physician indicates that the correction is good, it may not be necessary to force the heel down into the AFO, as this can lead to skin irritation, bracing discomfort, and reduced bracing compliance.
How long will each pair of AFOs last?	On average, an AFO will last 3-9 months. Some children grow rapidly and will need new AFO's within three months. As the child's growth rate slows, more time can be expected from each AFO.
How long will the bar last?	The bars are adjustable and should last approximately more than one pair of AFOs. The bar width should be periodically checked and adjusted to match the child's shoulder width or as recommended by your physician or orthotist.
Can I adjust the springs?	The springs are set, and should not be adjusted.
Can child walk in the <b>PONSETI® Abduction Bar</b> ?	It is not recommended that the child walk in the <b>Ponseti® Abduction Bar</b> . Walking with the bar attached puts the child at a high fall risk and may cause breakage to the brace which could negatively affect correction.
Can child walk in the <b>MITCHELL PONSETI® Move Bar</b> ?	It is not recommended that the child walk in the <b>Mitchell Ponseti® Move Bar</b> . Walking with the bar attached puts the child at a high fall risk and may cause breakage to the brace which could negatively affect correction.

## 8 GLOSSARY

Term	Meaning
Ankle Foot Orthosis (AFO)	An AFO is a brace, usually made of plastic, that is worn on the lower leg and foot to support the ankle, hold the foot and ankle in the correct position
Casting	A procedure to keep the clubfoot in place after manipulation. Padding is wrapped around the foot, then plaster or fiberglass is applied, and the cast is molded into place.
Clubfoot (talipes equinovarus)	Clubfoot (Talipes equinovarus) is a common birth defect in which the foot is twisted out of shape or position. In clubfoot, the foot appears twisted and can even look as if it's upside down. Despite its appearance, clubfoot itself doesn't cause any discomfort or pain. Treatment is usually successful and includes stretching and casting (Ponseti Method) or stretching and taping (French method). Sometimes, surgery is needed.
Dorsiflexion	Dorsiflexion is the flexion of the foot in an upward direction, from the ankle.
Equinus	Equinus is a condition in which the upward bending motion of the ankle joint is limited. Someone with equinus lacks the flexibility to bring the top of the foot toward the front of the leg. Equinus can occur in one or both feet.
Foot Abduction Brace	Abduction is a motion that pulls a structure or part away from the midline of the body. An abduction brace keeps the foot in proper alignment with the midline.
Orthotist	An orthotist is a healthcare professional who makes and fits braces and splints (orthoses) for people who need added support for body parts that have been weakened by injury, disease, or disorders of the nerves, muscles, or bones
Ponseti Method	The Ponseti Method is a non-surgical clubfoot treatment process developed by Dr. Ignacio Ponseti. The method utilizes gentle manipulation and casting to gradually move the foot into the correct position.
Pressure Saddles	A pressure saddle is a piece of soft, flexible rubber that attaches to the middle strap of the AFO to provide additional relief from the pressure of the strap.
Tenotomy	An Achilles tendon release, also known as tenotomy, is a surgical procedure that involves cutting through the Achilles tendon to allow for final dorsiflexion of the foot. The procedure is used to correct the equinus deformity seen in clubfoot.

## 9 SYMBOL LEGEND

Symbol	Meaning	Symbol	Meaning
	Ankle-Foot Orthotic (AFO)		Caution
	Authorized Representative in the European Community		Date of manufacture
	Manufacturer		Importer of record
	Single Patient, Multiple Use		Non-sterile
	Indicates Item is a Medical Device		Catalog number
	CE Mark (in compliance with Regulation (EU) 2017/745 of the European Parliament and of the council of 5 April 2017 on Medical Devices)		Batch code
	Information considered important, but not hazard-related		Electronic Instructions for Use

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To: MD ORTHOPAEDICS Date: 03/26/24 Proof # 2-PDF

Part No: IFU MD-124.3 REV 10 Env. ID: M8089 (WGS) Prod. # 04J-451-01

Reference: Users Instructions - Slavic, April 2023 PO # Artwork Ticket

Flat Size: N/A Finished Size: 7.00000" X 9.00000"

Ink Color(s): BLACK - 1/1

One-Sided /  Two-sided  Booklet Page Count: 136 Revision Date: 03/26/24

**1. Please proof the preceding pages: READ CAREFULLY & CHECK FOR ERRORS.**

*Mark any corrections directly on proof and return to Robinson Printing.*

We cannot be responsible for changes not made in writing. Please check for errors, omissions, etc. Signing below releases Robinson Printing from any liability regarding accuracy and/or correctness of design, grammar, spelling or information of this project in its entirety. Responding by email with approval to print also constitutes a binding contract.

**2. Please check appropriate box:**

- APPROVED TO PRINT -- NO CORRECTIONS
- ARTWORK APPROVAL ONLY -- DO NOT PRINT
- HARD COPY PROOF REQUESTED
- NEEDS CORRECTIONS -- requires new proof

**3. Please return this signed portion with changes, if any, so we may proceed with your job.**

**See Robinson Printing/Steven Label Terms and Conditions online at :**

**<https://stat.stevenlabel.com/termsandconditions.asp>**

**I take full responsibility for the approved artwork.**

Name (printed) \_\_\_\_\_

Signed \_\_\_\_\_ Date \_\_\_\_\_

**Comments:**

Artist: JSW CSR: KGO

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Copyright © 2023 MD Orthopaedics, Inc.